One-Step RNA Removal Kit

Introduction

One-Step RNA Removal Kit uses Ribonuclease A immobilized magnetic beads to efficiently remove RNA from the sample using a single step protocol. Ribonuclease A is an endoribonuclease that originates from the bovine pancreas. RNase A is a single chain polypeptide with a molecular mass of 13.7 kDa. RNase A is an endoribonuclease that unspecific ally degrades ribonucleic acid (RNA) into smaller components. The magnetic bead immobilized with RNase A can efficiently remove RNA from biological samples with no nucleases remaining in the solution due to the nuclease stably and covalently conjugated with the magnetic beads.


Applications

- For analytical purposes
- Plasmid and genomic DNA preparation
- · For cell cycle analysis by flow cytometry and propidium iodide (PI) staining
- Removal of RNA from recombinant protein preparations
- Ribonuclease protection assays. Used in conjunction with RNase T1
- Mapping single-base mutations in DNA or RNA

Features and Advantages

- Efficient one-tube and extraction-free protocol (Fig.1).
- Ultrafast: Process 96 samples in less than 30 minutes with <10-second Hands-on Time
- Nuclease Recovered at the end of the reaction thereby can be reused.
- Easy separation of the endonuclease from the reaction.
- · Stability of the immobilized nuclease increases.
- · Cost-effective: Eliminates columns, filters, laborious, organic reagents, and minimal plasticware required.
- High throughput: Compatible with many different automated liquid handling systems.

Workflow



Formulation: Liquid (Supplied in 50 mM Tris-HCl, pH 8.0, 50% Glycerol.)

Activity: 1 µl Magnetic Beads will digest 0.5µg of yeast RNA (Sigma, Catalog # R9001) in 50 mM Tris-HCl, pH 8.0 in 15 minutes at 37 °C. Shipping: Shipped at ambient temperature. Upon receipt, store nuclease magnetic Beads at -20°C. Aliquot to avoid repeated freezing and thawing.

Products	Catalog # AU-101	Catalog # AU-102
BcMag [™] One-Step RNA Removal Kit	1 ml	3 ml

PROTOCOL

A. Accessory equipment

Magnetic Rack

Item	Source
Magnetic Rack for centrifuge tube	BcMag Rack-2 for holding two individual 1.5 ml centrifuge tubes
	(Bioclone, Cat. # MS-01)



Magnetic Beads Make Things Simple

Instruction Manual

2

** Based on sample volume, the user can choose one of the following magnetic Racks	BcMag Rack-6 for holding six individual 1.5 ml centrifuge tubes (Bioclone, Cat. # MS-02) BcMag Rack-24 for holding twenty-four individual 1.5-2.0 ml centrifuge tubes (Bioclone, Cat. # MS-03) BcMag Rack-50 for holding one 50 ml centrifuge tube, one 15 ml centrifuge tube, and four individual 1.5 ml centrifuge tubes (Bioclone, Cat. # MS-04)
BcMag 96-well Plate Magnetic Rack.	BcMag 96-well Plate Magnetic Rack (side-pull) compatible with 96-well PCR plate and 96-well microplate or other compatible Racks (Blioclone, Cat#: MS-05)

B. Procedure

- Do not use buffers containing organic solvents.
- Typically, the bead is added directly into any standard buffer at the desired amount of the beads based on the concentration of the RNA (1 µl Magnetic Beads will digest 10µg of yeast RNA).

Table 1. Working Condition

Conditions	Optimal Function	Functional Range	Inhibitory Action
Temperature	60 °C	15-70 °C	N/A
pН	7.6	6.0-10.0	N/A
Zn^{2+}			Yes
Cu ²⁺	N/A	N/A	Yes
0-100 mM NaCl	Preferentially cleave single-stranded and double-stranded RNA and the RNA strand in RNA-DNA hybrids.	N/A	
>300 mM NaCl	Preferentially cleave single-stranded RNA		
Guanidine HCl	N/A	N/A	4M + 0.1 M 2-mercaptoethanol
SDS			Yes

1. Shake the bottle to resuspend the Magnetic beads until it is homogeneous entirely.

IMPORTANT! It is essential to mix the beads before dispensing. Do not allow the beads to sit for more than 2 minutes before dispensing. Resuspend the magnetic beads every 2 minutes.

- 2. Add an appropriate amount of the magnetic Beads to a reaction.
- 3. Mix the sample with beads for 1-2 minutes by slowly pipetting up and down 20-25 times or Vortex the sample for 2 minutes at 2000 rpm.
- 4. Incubate at 37°C with continuous rotation for 15 minutes.
- 5. Place the sample plate or tube on the magnetic Rack for 30 seconds or until the solution is clear. (Option: centrifuge at 3500 rpm for 45 seconds)
- 6. Transfer the supernatant to a clean plate /tube while the sample plate remains on the magnetic separation plate. The sample is ready for downstream applications.

Reference

1. Molecular Cloning, A Laboratory Manual (3rd ed). Cold Spring Harbor Laboratory Press (Cold Spring Harbor, NY), Volume 1, 1.78-1.79 (2001).

Related products

Products and Catalog Number				
Genomic DNA and RNA Purification				
One-Step Mammalian Cell DNA Purification Kit, Cat. No. AA101	One-Step Saliva Viral RNA-DNA Purification Kit, Cat. No. AR101			
Cell-Free DNA Purification Kit, Cat. No AC101	Bone-Teeth DNA Purification Kit, Cat. No. AB101			
One-Step FFPE & FNA DNA purification Kit, Cat. No. AJ-101	Rootless Hair DNA Purification Kit, Cat. No. AD101			
One-Step Bacteria DNA Purification Kit, Cat. No. AE101	One-Step Buccal Cell DNA Purification Kit, Cat. No. AG101			
One-Step Blood DNA Purification Kit, Cat. No. AF101	One-Step Touch DNA Purification Kit, Cat. No. AS101			
One-Step Fungi & Yeast DNA Purification Kit, Cat. No. AL101	Sexual Assault Casework DNA Purification Kit, Cat. No. AT101			
One-Step Insect DNA Purification Kit, Cat. No. AM101	One-Step Fingerprint DNA Purification Kit, Cat. No. AZ101			
One-Step Mouse Tail DNA Purification Kit, Cat. No. AN101	One-Step Dandruff DNA Purification Kit, Cat. No. AAA101			



Magnetic Beads Make Things Simple

Instruction Manual

3

One-Step Plant DNA Purification Kit, Cat. No. AQ101	Quick mRNA Purification Kit, Cat. No. MMS101		
DNA & RNA Sample Preparation			
One-Step NGS Cleanup Kit, Cat. No. AO101	One-Step DNA-RNA Removal Kit, Cat. No. CA103		
One-Step RNA Removal Kit, Cat. No. AU101	One-Step DNA/RNA Cleanup Kit, Cat. No. AH101		
One-Step PCR Cleanup Kit, Cat. No. AP101	One-Step Sequencing Cleanup Kit, Cat. No. AI101		
Quick Oligo-DNA Conjugation Kit, Cat. No. CA101	One-Step Fluorescent Labeling Cleanup Kit, Cat. No. AK101		
One-Step DNA-RNA Removal Kit, Cat. No. AV101	One-Step Single-Stranded DNA Removal Kit, Cat. No. AW101		
One-Step PCR Inhibitor Removal Kit, Cat. No. AX101	Pure Miniprep Plasmid DNA Purification Kit, Cat. No. AY101		

Phone: +1-858-909-0079 **Bioclone Inc.** Fax: +1-858-909-0057